

MR Neuro OR Tabletop for IMRIS Surgical Theatre

Specialty OR Tabletop for Neuro Surgery on TruSystem™ 7500 OR Table



TruSystem™ 7500 MR Neuro OR Table with floor mounted column and MR conditional* OR Tabletop

Enhanced Workflow for IMRIS Surgical Theatre

Perform safe and complex neuro surgery with Trumpf Medical's MR Neuro OR Tabletop.

- Optimal surgical site access due to individual two-segmented tabletop.
- Anatomical adjusted head positioning with special head fixation devices.
- Ideal intraoperative imaging access due to lean design.
- Enhanced flexibility with inter-changeable TruSystem™ 7500 OR Table System.

The MR Neuro OR Tabletop supports the OR team for efficient clinical outcomes.

*MR conditions are described in the MR Neuro Tabletop IFU.

Precise and Flexible



Optimal Intraoperative Access

Adaptable Head-Fixation Devices (HFD) 100/200 by IMRIS for flexible head positioning.



Enhanced Workflow

TruSystem™ 7500 MR Neuro OR Table and IMRIS VISIUS iMRI System are ideally harmonized in scanning.



Advanced Table Configuration

Easy attachable side rails and components for positioning the table precisely into the iMRI scanner.



General to Specialty

Transform your OR to accommodate multiple specialties with one single column for the TruSystem™ 7500 OR Table and a wide range of tabletops.

Technical Details

Material No.	MR Neuro OR Tabletop: EU # 2064887 US # 2064888 JP # 2064889 OR Column TruSystem 7500 MR Neuro OR Table: # 2064886
Dimensions of Tabletop (L x W)	2412 mm x 520 mm 95" x 20.5"
Side Rail Dimensions	EU: 25 mm x 10 mm 1" x 0.4" US: 28.6 mm x 9.5 mm 1.1" x 0.4" JP: 32 mm x 9 mm 1.3" x 0.4"
Longitudinal Slide	Stationary Column: 472 mm 18.6"
Net Weight	ca. 197 kg 434 lbs
Maximum Load	Patient Weight: 180 kg 397 lbs Fixed Side Rail: 50 Nm Clipped Side Rail: 25 Nm

Not all products/ options are available in all countries. For further information about our products or services, please contact your local Trumpf Medical representative or visit our webpage:

www.trumpfmedical.com

Trumpf Medical reserves the right to make changes without notice in design, specifications and models. The only warranty Trumpf Medical makes is the express written warranty extended on the sale or rental of its products. ©2018 TRUMPF Medizin Systeme GmbH + Co. KG, ALL RIGHTS RESERVED.

TRUMPF Medizin Systeme GmbH + Co. KG
Benzstraße 26 · 82178 Puchheim · Telefon +49 (0) 89 80907-0 · E-Mail info@trumpfmedical.com

2069302, 1 February 2018